

Anti-Wnt2b Rabbit antibody



Catalog Number: 168894

Product name

Anti-Wnt2b Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Wnt2b Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human Wnt2b

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC


Dilutions

1:1000

1:20

1:20

Validations



Western blot detection of Wnt2b in K562, C6, 3T3, HeLa cell lysates using Wnt2b Rabbit mAb(1:1000 diluted). Predicted band size: 44kDa. Observed band size: 44kDa.