

# Human PPID/CYP40/CYPD (His tag) recombinant protein



Catalog Number: 519894

## General Information

### Gene Name Synonym

Peptidyl-Prolyl Cis-Trans Isomerase D, PPIase D, 40 kDa Peptidyl-Prolyl Cis-Trans Isomerase, Cyclophilin-40, CYP-40, Cyclophilin-Related Protein, Rotamase D, PPID, CYP40, CYPD

### Protein Construction

A DNA sequence encoding the Human PPID/CYP40/CYPD Met1-Ala370 is expressed with a 6His tag at the N-terminus, 6His tag at the C-terminus.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu$ g protein as determined by the LAL method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2  $\mu$ m filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human PPID/CYP40/CYPD (His tag) recombinant protein