

GSTT2B antibody



Catalog Number: 111194

Product name

GSTT2B antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

GSTT2B Rabbit Polyclonal antibody. Positive IF detected in HepG2. Positive WB detected in HepG2 cells, human brain tissue, human liver tissue. Observed molecular weight by Western-blot: 28kd

Preparation

This antibody was obtained by immunization of GSTT2B recombinant protein (Accession Number: NM_001080843). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:111194(GSTT2B antibody) at dilution of 1:500



Immunofluorescent analysis of HepG2 cells, using GSTT2B antibody Catalog No:111194 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).