

LENG8 antibody



Catalog Number: 112194

Product name

LENG8 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

LENG8 Rabbit Polyclonal antibody. Positive WB detected in A2780 cells. Positive FC detected in HEK-293 cells. Observed molecular weight by Western-blot: 88 kDa

Preparation

This antibody was obtained by immunization of LENG8 recombinant protein (Accession Number: NM_052925). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, FC

Dilutions


Recommended Dilution:

WB: 1:500-1:5000

Validations



A2780 cells were subjected to SDS PAGE followed by western blot with Catalog No:112194(LENG8 antibody) at dilution of 1:1200



1X10⁶ HEK-293 cells were stained with 0.2ug LENG8 antibody (Catalog No:112194, red) and control antibody (blue). Fixed with 4% PFA blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1500.