

Anti-Claudin 3 Rabbit antibody



Catalog Number: 167194

Product name

Anti-Claudin 3 Rabbit antibody

Specificity

Human

Antibody description

Claudin 3 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human Claudin 3

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC, IP

Dilutions

WB: 1/1000


IHC: 1/50

ICC/IF: 1/20

FC: 1/20

IP: 1/20

Validations



Western blot detection of Claudin 3 in MCF7 cell lysates using Claudin 3 Rabbit pAb(1:1000 diluted). Predicted band size:23KDa. Observed band size:20KDa.