

# Anti-PKM2 (1A7) Mouse antibody



Catalog Number: 168194

## Product name

Anti-PKM2 (1A7) Mouse antibody

## Specificity

Human, Mouse, Rat, Monkey

## Antibody description

PKM2 (1A7) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human PKM2 protein fragments expressed in E.coli.

## Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


WB, ICC/IF

## Dilutions


WB: 1/1000

;ICC: 1/400

## Validations



Immunocytochemistry staining of HeLa cells fixed with 4% Paraformaldehyde and using anti-PKM2 mouse mAb (dilution 1:400).



Western blot detection of PKM2 in COS7, PC12, C6, 3T3 and HeLa cell lysates using PKM2 mouse mAb (1:1000 diluted). Predicted band size: 60KDa. Observed band size: 60KDa.