

# Icln antibody

Catalog Number: 111594

## Product name

Icln antibody

## Immunogen

[Human CLNS1A Recombinant protein \(GST tag\)](#)

## Specificity

Human, Dog; other species not tested.

## Antibody description

Icln Rabbit Polyclonal antibody. Positive IF detected in HEK-293 cells. Positive WB detected in HEK-293 cells, MDCK cells. Observed molecular weight by Western-blot: 36-42 kDa

## Preparation

This antibody was obtained by immunization of Icln recombinant protein (Accession Number: NM\_001311199). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:111594(CLNS1A Antibody) at dilution of 1:1000



Immunofluorescent analysis of (-20oc Ethanol) fixed HEK-293 cells using Catalog No:111594(CLNS1A Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)