

MED8 antibody



Catalog Number: 112594

Product name

Polyclonal

MED8 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse; other species not tested.

Applications

Antibody description

ELISA, WB

MED8 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, mouse liver tissue, SGC-7901 cells. Observed molecular weight by Western-blot: 29 kDa

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Preparation

This antibody was obtained by immunization of MED8 recombinant protein (Accession Number: BC010543). Purification method: Antigen affinity purified.

Validations

35kd→
25kd→
18kd→
14kd→



Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:112594(MED8 antibody) at dilution of 1:1000

Clonality