

BNIP3 Polyclonal Antibody



Catalog Number: 164594

Product name

BNIP3 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to BNIP3

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-164 of human BNIP3 (NP_004043.3).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:200

Validations



Western blot - BNIP3 Polyclonal Antibody


Western blot analysis of extracts of various cell lines, using BNIP3 Antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL



Basic Kit .

Immunohistochemistry - BNIP3 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse kidney using BNIP3 antibody at dilution of



1:100 (40x lens).


BNIP3 Polyclonal Antibody



Catalog Number: 164594

Immunohistochemistry - BNIP3 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat heart using BNIP3 antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - BNIP3 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human liver injury using BNIP3 antibody at dilution of 1:100 (40x lens).



Immunofluorescence - BNIP3 Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using BNIP3 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.