

# USP18 antibody

Catalog Number: 116594

## Product name

USP18 antibody

## Specificity

Human; other species not tested.

## Antibody description

USP18 Rabbit Polyclonal antibody. Positive WB detected in MCF-7. Observed molecular weight by Western-blot: 43 kDa, 70 kDa

## Preparation

This antibody was obtained by immunization of USP18 recombinant protein (Accession Number: NM\_017414). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: N/A

## Validations



WB result of USP18 antibody(Catalog No:116594) with MCF-7 cells (Control and CRISPR USP18).