

CSF2RA/CD116 antibody



Catalog Number: 109594

Product name

CSF2RA/CD116 antibody

Specificity

Human; other species not tested.

Antibody description

CSF2RA/CD116 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, human placenta tissue. Positive IHC detected in human placenta tissue. Observed molecular weight by Western-blot: 48kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: XM_011545619). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:109594(CSF2RA antibody) at dilution of 1:600



Immunohistochemical of paraffin-embedded human placenta using Catalog No:109594(CSF2RA antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human placenta using Catalog No:109594(CSF2RA antibody) at dilution of 1:50 (under 40x lens)