

# Anti-ICAM2 Rabbit antibody



Catalog Number: 167994

## Product name

Anti-ICAM2 Rabbit antibody

## Specificity

Mouse

## Antibody description

ICAM2 Rabbit monoclonal antibody

## Preparation

Antigen: Recombinant protein of mouse ICAM2

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, IP


## Dilutions

WB: 1/1000

IHC: 1/500

IP: 1/20

## Validations



Western blot detection of ICAM2 in HeLa,A549,HL-60,U2OS cell lysates using ICAM2 Rabbit mAb(1:1000 diluted).Predicted band size:31kDa.Observed band size:50-60kDa.