

Anti-Asparagine Synthetase Rabbit antibody



Catalog Number: 168994

Product name

Anti-Asparagine Synthetase Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Human, Mouse, Rat

Applications

WB, IHC-P, ICC/IF

Antibody description

Asparagine Synthetase Rabbit monoclonal antibody

Dilutions

1:1000

1:20

1:20-1:50


Preparation

Antigen: A synthetic peptide of human Asparagine synthetase

Validations

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.



Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Western blot detection of Asparagine Synthetase in K562, HeLa cell lysates using Asparagine Synthetase Rabbit mAb(1:1000 diluted). Predicted band size:64kDa. Observed band size:64kDa.