

SLAM/CD150 antibody



Catalog Number: 115294

Product name

SLAM/CD150 antibody

Specificity

Human; other species not tested.

Antibody description

SLAM/CD150 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, human placenta tissue. Observed molecular weight by Western-blot: 70 kDa

Preparation

This antibody was obtained by immunization of SLAM/CD150 recombinant protein (Accession Number: BC012602). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:115294(SLAMF1 Antibody) at dilution of 1:300