

Rabbit anti Myc-Tag pAb



Catalog Number: 166294

Product name

Rabbit anti Myc-Tag pAb

Antibody description

Polyclonal antibody to Rabbit anti Myc-Tag

Preparation

Antigen: Recombinant protein of MYC-tag

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:1000 - 1:10000

Validations



Western blot - Rabbit anti Myc-Tag pAb

Western blot analysis of over-expressed Myc-tagged protein in 293T cell using Myc-tag antibody at different dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.