

Gelsolin antibody



Catalog Number: 107294

Product name

Gelsolin antibody

Recommended Dilution:

WB: 1:1000-1:10000

Immunogen

IHC: 1:20-1:200

[Human Gelsolin Recombinant protein \(His tag\)](#)


Validations

Specificity

Human; other species not tested.

Antibody description

Gelsolin Mouse Monoclonal antibody. Positive IHC detected in human colon tissue, human colon cancer tissue. Positive WB detected in human plasma tissue, A549 cells. Observed molecular weight by Western-blot: 80-83 kDa



Preparation

This antibody was obtained by immunization of Gelsolin recombinant protein (Accession Number: NM_000177). Purification method: Thiophilic affinity chromatograph.

human plasma tissue were subjected to SDS PAGE followed by western blot with Catalog No:107294(Gelsolin Antibody) at dilution of 1:2000

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal


Ig Type

Mouse IgA

Applications

ELISA, WB, IHC


Dilutions



Immunohistochemistry of paraffin-embedded human colon tissue slide using Catalog No:107294(Gelsolin Antibody) at dilution of 1:200 (under 10x lens).

Gelsolin antibody

Catalog Number: 107294



Immunohistochemistry of paraffin-embedded human colon tissue slide using Catalog No:107294(Gelsolin Antibody) at dilution of 1:200 (under 40x lens).