

FOXD4 antibody



Catalog Number: 110694

Product name

Polyclonal

FOXD4 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

Applications

Antibody description

ELISA, WB

FOXD4 Rabbit Polyclonal antibody. Positive WB detected in NIH/3T3 cells, Jurkat cells, RAW 264.7 cells, U-937 cells. Observed molecular weight by Western-blot: 65-70 kDa

Dilutions


Recommended Dilution:

WB: 1:500-1:5000

Preparation

This antibody was obtained by immunization of FOXD4 recombinant protein (Accession Number: NM_207305). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

NIH/3T3 cells were subjected to SDS PAGE followed by western blot with Catalog

No:110694(FOXD4 Antibody) at dilution of 1:1000

Clonality