

Anti-NANOG antibody



Catalog Number: 141694

Product name

Anti-NANOG antibody

Specificity

Human NANOG / Nanog homeobox

Antibody description

Rabbit Polyclonal to Human NANOG

Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the C-terminus of the Human NANOG / Nanog homeobox, and purified by antigen affinity chromatography.

Formulation

0.2 μ m filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IHC-P

Dilutions

IHC-P: 0.1-2 μ g/mL

Validations



NANOG / Nanog homeobox Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry