

DLGAP5 Polyclonal Antibody



Catalog Number: 161694

Product name

DLGAP5 Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to DLGAP5

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 547-846 of human DLGAP5 (NP_055565.3).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000


IF 1:50 - 1:200

Validations



Western blot - DLGAP5 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using DLGAP5 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Mouse IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL



Basic Kit .Exposure time: 90s.

Immunofluorescence - DLGAP5 Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using DLGAP5 antibody . Blue: DAPI for nuclear staining.