

Anti-SLC38A2 antibody

Catalog Number: 175694

Product name

Anti-SLC38A2 antibody

Specificity

Human, Mouse, Rat, Dog, Pig, Cow, Horse, Sheep

Antibody description

Rabbit polyclonal antibody to SLC38A2

Preparation

This antigen of this antibody was klh conjugated synthetic peptide derived from human slc38a2/snat2 21-150/506

Formulation

Liquid, 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

Storage

Store at -20°C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4°C.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB:1:500-2000

IHC-P:1:400-800

Validations




Fig1: Sample:; Heart (Mouse) Lysate at 40 ug; Primary: Anti-SLC38A2 at 1/300 dilution; Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution; Predicted band size: 56 kD; Observed band size: 56 kD




Fig2: Sample:; Intestine (Mouse) Lysate at 40 ug; Primary: Anti-SLC38A2 at 1/300 dilution; Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution; Predicted band size: 56 kD; Observed band size: 56 kD




Fig3: Paraformaldehyde-fixed, paraffin embedded

Anti-SLC38A2 antibody



Catalog Number: 175694

(rat brain); Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15min; Block endogenous peroxidase by 3% hydrogen peroxide for 20 minutes; Blocking buffer (normal goat serum) at 37°C for 30min; Antibody incubation

with (SLC38A2) Polyclonal Antibody, Unconjugated at 1:400 overnight at 4°C, followed by a conjugated secondary (sp-0023) for 20 minutes and DAB staining.