

Anti-Phospho-p44/42 MAPK (ERK1/2) (Thr202/Tyr204) Rabbit antibody



Catalog Number: 166694

Product name

Anti-Phospho-p44/42 MAPK (ERK1/2)
(Thr202/Tyr204) Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Phospho-p44/42 MAPK (ERK1/2) (Thr202/Tyr204)
Rabbit polyclonal antibody

Preparation

Antigen: Synthetic peptide corresponding to
human Phospho-p44/42 MAPK (Erk1/2)
(Thr202/Tyr204) protein

Formulation

PBS(pH 7.4) containing with 0.03% Proclin300
and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 1/1000

Validations



Western blot analysis of Phospho-p44/42 MAPK (Erk1/2) (Thr202/Tyr204) from EGF-treated A549 cells and Hydroxyurea-treated NIH/3T3 cells, using Phospho-p44/42 MAPK (Erk1/2) (Thr202/Tyr204) Rabbit pAb (166694, 1:1000 diluted, upper) and p44/42 MAPK (Erk1/2) Mouse mAb (201246-4F3, 1:1000 diluted, lower).