

HEBP1 antibody



Catalog Number: 111394

Product name

HEBP1 antibody

ELISA, WB, IHC

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

IHC: 1:20-1:200

Specificity


Human, Mouse, Rat; other species not tested.

Antibody description

HEBP1 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, HeLa cells, HepG2 cells, human brain tissue, rat liver tissue. Positive IHC detected in human liver cancer tissue.

Observed molecular weight by Western-blot: 21 kDa

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:111394(HEBP1 antibody) at dilution of 1:500

Preparation

This antibody was obtained by immunization of HEBP1 recombinant protein (Accession Number: NM_015987). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT


Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using Catalog No:111394(HEBP1 Antibody) at dilution of 1:50 (under 40x lens)