

Lumican antibody

Catalog Number: 112394

Product name

Lumican antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Lumican Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 65-90 kDa

Preparation

This antibody was obtained by immunization of Lumican recombinant protein (Accession Number: NM_002345). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:112394(LUM antibody) at dilution of 1:800