

Anti-ZNF75D antibody



Catalog Number: 176394

Product name

Anti-ZNF75D antibody

Ig Type

Rabbit IgG

Specificity

Human

Applications

WB, IHC-P

Antibody description

Rabbit polyclonal antibody to ZNF75D

Dilutions

WB: 1:500-1:2000

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human znf75. aa range:296-345

IHC-P: 1:100-1:300

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Validations




Fig1:: Western Blot analysis of various cells using ZNF75 Polyclonal Antibody diluted at 1: 500. Secondary antibody (catalog#: HA1001) was diluted at 1:20000