

Anti-ZNF75D antibody



Catalog Number: 176394

Product name

Anti-ZNF75D antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to ZNF75D

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human znf75. aa range:296-345

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations



Fig1.; Western Blot analysis of various cells using ZNF75 Polyclonal Antibody diluted at 1: 500. Secondary antibody (catalog#: HA1001) was diluted at 1:20000