

# B4GALT3 antibody



Catalog Number: 108394

---

## Product name

B4GALT3 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

B4GALT3 Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer.

## Preparation

This antibody was obtained by immunization of B4GALT3 recombinant protein (Accession Number: XM\_024450540). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC

## Dilutions

Recommended Dilution:

IHC: 1:20-1:200

## Validations



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:108394(B4GALT3 antibody) at dilution of 1:50 (under 10x lens)