

# ADIPOR2 antibody

Catalog Number: 107794

## Product name

ADIPOR2 antibody

## Immunogen

[Human ADIPOR2 Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ADIPOR2 Rabbit Polyclonal antibody. Positive WB detected in L02 cells, HeLa cells, mouse uterus tissue. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 44-50 kDa

## Preparation

This antibody was obtained by immunization of ADIPOR2 recombinant protein (Accession Number: XM\_047429546). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:107794(ADIPOR2 antibody) at dilution of 1:500



IP Result of anti-ADIPOR2 (IP:Catalog No:107794, 4ug; Detection:Catalog No:107794 1:500) with HeLa cells lysate 2800ug.