

Anti-Csf2ra antibody



Catalog Number: 175005

Product name

Anti-Csf2ra antibody

Specificity

Mouse

Antibody description

Rabbit polyclonal antibody to Csf2ra

Preparation

This antigen of this antibody was synthetic peptide within n-terminal mouse gm-csf-r-alpha (cd116).

Formulation

Liquid, 1*PBS (pH7.4), 0.2% BSA, 40% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ICC, IHC-P, FC

Dilutions

WB: 1:500-1:1000

ICC: 1:100-1:200

FC: 1:10-1:100

Validations




Fig1: Western blot analysis on F9 cell lysates using anti- GM-CSF-R-alpha polyclonal antibody.