

SCN10A antibody



Catalog Number: 115005

Product name

SCN10A antibody

Specificity

Human; other species not tested.

Antibody description

SCN10A Rabbit Polyclonal antibody. Positive IF detected in SH-SY5Y cells.

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_006514).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol
pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IF

Dilutions

Recommended Dilution:

IF: 1:10-1:100

Validations



Immunofluorescent analysis of SH-SY5Y cells using Catalog No:115005(SCN10A Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)