

# ZNF571 antibody

Catalog Number: 117005

## Product name

ZNF571 antibody

## Immunogen

[Human ZNF571 Recombinant protein \(His tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

ZNF571 Rabbit Polyclonal antibody. Positive IHC detected in human testis tissue. Positive IP detected in MCF-7 cells. Positive WB detected in MCF-7 cells, HeLa cells. Observed molecular weight by Western-blot: 65 kDa

## Preparation

This antibody was obtained by immunization of ZNF571 recombinant protein (Accession Number: NM\_001290314). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC, IP

## Dilutions


Recommended Dilution:

WB: 1:200-1:1000


IP: 1:200-1:1000

IHC: 1:20-1:200

## Validations



MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:117005(ZNF571 Antibody) at dilution of 1:300




IP Result of anti-ZNF571 (IP:Catalog No:117005, 4ug; Detection:Catalog No:117005 1:300) with MCF-7 cells lysate 1200ug.

# ZNF571 antibody

Catalog Number: 117005

---



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:117005(ZNF571 Antibody) at dilution of 1:50 (under 40x lens)