

ERC1 antibody



Catalog Number: 110405

Product name

ERC1 antibody

Specificity

Human; other species not tested.

Antibody description

ERC1 Rabbit Polyclonal antibody. Positive IF detected in MCF-7 cells. Positive WB detected in HeLa cells, COLO 320 cells, human testis tissue, Jurkat cells. Observed molecular weight by Western-blot: 124-135 kDa

Preparation

This antibody was obtained by immunization of ERC1 recombinant protein (Accession Number: XM_011520941). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110405(ERC1 antibody) at dilution of 1:1000



Immunofluorescent analysis of MCF-7 cells, using ERC1 antibody Catalog No:110405 at 1:50 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).