

Anti-IL20RA/IL-20RA antibody



Catalog Number: 103805

Product name

Anti-IL20RA/IL-20RA antibody

Immunogen

[Human IL20RA/IL-20RA \(His Tag\) recombinant protein](#)

Specificity

Human IL20RA / IL-20RA

No cross-reactivity in ELISA with Human IFN γ R1; Human IL10RB

Antibody description

Rabbit monoclonal to IL20RA/IL-20RA

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human IL20RA / IL-20RA (rh IL20RA / IL-20RA; NP_055247.21; Met 1-Lys 250).

Formulation

0.2 μ m filtered solution in PBS with 5% trehalose

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC-P


Dilutions

ELISA: 0.5-1 μ g/mL

This antibody can be used at 0.5-1 μ g/mL with the appropriate secondary reagents to detect Human IL20Ra.

IHC-P: 1-10 μ g/mL

Validations



IL20RA / IL-20RA Antibody, Rabbit MAb, Immunohistochemistry

Immunochemical staining of human IL20Ra in human skin with rabbit monoclonal antibody (5 μ g/mL, formalin-fixed paraffin embedded sections). Positive staining was localized to stratified squamous epithelium.