

VTA1 antibody

Catalog Number: 116805

Product name

VTA1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

VTA1 Rabbit Polyclonal antibody. Positive IF detected in HeLa cells. Positive WB detected in human kidney tissue, human heart tissue. Observed molecular weight by Western-blot: 45kd

Preparation

This antibody was obtained by immunization of VTA1 recombinant protein (Accession Number: NM_016485). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IF: 1:10-1:100

Validations



human kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:116805(VTA1 antibody) at dilution of 1:500



Immunofluorescent analysis of HeLa cells, using VTA1 antibody Catalog No:116805 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).