

DNAJC15 antibody



Catalog Number: 110105

Product name

DNAJC15 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

DNAJC15 Rabbit Polyclonal antibody. Positive IHC detected in human liver. Positive WB detected in mouse liver tissue, human brain tissue. Observed molecular weight by Western-blot: 18-20kd

Preparation

This antibody was obtained by immunization of DNAJC15 recombinant protein (Accession Number: NM_013238). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200


Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:110105(DNAJC15 antibody) at dilution of



1:1000
Immunohistochemical of paraffin-embedded human liver using Catalog No:110105(DNAJC15 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human liver using Catalog No:110105(DNAJC15 antibody) at dilution of 1:100 (under 40x lens)