

Anti-RSK1 p90 (6B9) Mouse antibody



Catalog Number: 168105

Product name

WB, IP

Anti-RSK1 p90 (6B9) Mouse antibody

Dilutions

Specificity

WB: 1/1000

Human, Mouse, Monkey

Antibody description

RSK1 p90 (6B9) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human p90RSK protein fragments expressed in E.coli.

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality


Monoclonal

Ig Type


Mouse IgG1

Applications

Validations



Immunoprecipitation analysis of HeLa cell lysates using p90RSK mouse mAb.



Western blot detection of p90RSK in COS7,3T3,Jurkat,HeLa and A549 cell lysates using p90RSK mouse mAb(dilution 1:1000).Predicted band size:90KDa.Observed band size:90KDa.