

# Anti-RSK1 p90 (6B9) Mouse antibody



Catalog Number: 168105

## Product name

Anti-RSK1 p90 (6B9) Mouse antibody

## Specificity

Human, Mouse, Monkey

## Antibody description

RSK1 p90 (6B9) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human p90RSK protein fragments expressed in E.coli.

## Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1


## Applications

WB, IP


## Dilutions

WB: 1/1000

## Validations



Immunoprecipitation analysis of HeLa cell lysates using p90RSK mouse mAb.



Western blot detection of p90RSK in COS7,3T3,Jurkat,HeLa and A549 cell lysates using p90RSK mouse mAb(dilution 1:1000).Predicted band size:90KDa.Observed band size:90KDa.