

Anti-COL13A1 antibody

Catalog Number: 176505

Product name

Anti-COL13A1 antibody

Specificity

Human, Mouse

Antibody description

Rabbit polyclonal antibody to COL13A1

Preparation

This antigen of this antibody was synthesized peptide derived from human col13a1. at aa range: 41-90

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1:; Western Blot analysis of 1,Hela 2,Mouse-kidney cells using primary antibody diluted at 1:2000(4°C overnight). Secondary antibody: Goat Anti-rabbit IgG IRDye 800(diluted at 1:5000, 25°C, 1 hour)