

Anti-Olfactory Marker Protein/OMP Protein Rabbit antibody



Catalog Number: 167505

Product name

Anti-Olfactory Marker Protein/OMP Protein Rabbit antibody

Specificity

Mouse, Rat

Antibody description

Olfactory Marker Protein/OMP Protein Rabbit polyclonal antibody

Preparation

Antigen: Recombinant protein of mouse Olfactory Marker Protein

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, IP


Dilutions

WB: 1/1000

IHC: 1/200

IP: 1/20

Validations



Western blot detection of Olfactory Marker Protein in Rat Brain, MCF7 cell lysates using Olfactory Marker Protein Rabbit pAb(1:1000 diluted). Predicted band size: 19KDa. Observed band size: 19KDa.