

# BIRC8 antibody



Catalog Number: 108505

## Product name

BIRC8 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

BIRC8 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, human brain tissue. Observed molecular weight by Western-blot: 30-35 kDa

## Preparation

This antibody was obtained by immunization of BIRC8 recombinant protein (Accession Number: NM\_033341). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions


Recommended Dilution:

WB: 1:200-1:2000

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:117150(BIRC8 antibody) at dilution of 1:200



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:117150(BIRC8 antibody) at dilution of 1:300