

Anti-Pellino 1 Rabbit antibody



Catalog Number: 168905

Product name

Anti-Pellino 1 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Pellino 1 Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human Pellino 1

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, IP


Dilutions

WB: 1/1000

ICC/IF: 1/20

IP: 1/20

Validations



Western blot detection of Pellino 1 in Rat Brain, Jurkat cell lysates using Pellino 1 Rabbit mAb(1:500 diluted). Predicted band size: 46kDa. Observed band size: 46kDa.