

CSRNP1 Polyclonal Antibody



Catalog Number: 165205

Product name

CSRNP1 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to CSRNP1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human CSRNP1 (NP_149016.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:100

Validations



Immunohistochemistry - CSRNP1 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat kidney using CSRNP1 antibody at dilution of 1:200



(40x lens).

Immunohistochemistry - CSRNP1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse brain using CSRNP1 antibody at dilution of 1:200 (40x lens).



Immunofluorescence - CSRNP1 Polyclonal Antibody

CSRNP1 Polyclonal Antibody



Catalog Number: 165205

Immunofluorescence analysis of A549 cells using CSRNP1 antibody .