

SERPINB8 antibody



Catalog Number: 115205

Product name

SERPINB8 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

SERPINB8 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, mouse lung tissue. Observed molecular weight by Western-blot: 45-48kd

Preparation

This antibody was obtained by immunization of SERPINB8 recombinant protein (Accession Number: NM_001031848). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:115205(SERPINB8 antibody) at dilution of 1:200