

Fascin antibody

Catalog Number: 107205

Product name

Fascin antibody

Specificity

Human; other species not tested.

Antibody description

Fascin Mouse Monoclonal antibody. Positive IHC detected in human tonsillitis tissue. Positive WB detected in SH-SY5Y cells, HeLa cells, K-562 cells. Observed molecular weight by Western-blot: 55 kDa

Preparation

This antibody was obtained by immunization of Fascin recombinant protein (Accession Number: NM_003088). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2a

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200


Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:107205(Fascin Antibody) at dilution of 1:2000



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:107205(Fascin Antibody) at dilution of 1:600 (under 10x lens). Heat mediated antigen retrieved with Tris-EDTA buffer(pH9).



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog

Fascin antibody



Catalog Number: 107205

No:107205(Fascin Antibody) at dilution of 1:600

(under 40x lens). Heat mediated antigen retrieved with Tris-EDTA buffer(pH9).