

ART5 antibody

Catalog Number: 108205

Product name

ART5 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

ART5 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue. Positive IHC detected in human heart. Observed molecular weight by Western-blot: 35-40kd

Preparation

This antibody was obtained by immunization of ART5 recombinant protein (Accession Number: NM_053017). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:108205(ART5 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human heart using Catalog No:108205(ART5 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human heart using Catalog No:108205(ART5 antibody) at dilution of 1:50 (under 40x lens)