

KIAA0319 Polyclonal Antibody



Catalog Number: 162605

Product name

KIAA0319 Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to KIAA0319

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 21-300 of human KIAA0319 (NP_055624.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:200 - 1:2000

Validations



Western blot - KIAA0319 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using KIAA0319 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 5s.