

# GATA3 Monoclonal Antibody



Catalog Number: 164605

## Product name

GATA3 Monoclonal Antibody

## Specificity

Human

## Antibody description

Monoclonal antibody to GATA3

## Preparation

Antigen: A synthetic peptide of human GATA3

## Formulation

PBS with 0.02% sodium azide, pH7.3.

## Storage

Store at 4°C. Avoid freeze / thaw cycles.

## Clonality

Monoclonal

## Ig Type

Mouse IgG

## Applications

WB, IF, FC


## Dilutions

WB 1:500 - 1:2000

IF 1:20 - 1:50

FC 1:20 - 1:50

## Validations



## Western blot - GATA3 Monoclonal Antibody

Western blot analysis of extracts of various cell lines, using GATA3 antibody. Secondary antibody: HRP Goat Anti-Mouse IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.