

# Anti-Collagen I Rabbit antibody



Catalog Number: 167705

## Product name

Anti-Collagen I Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Collagen I Rabbit monoclonal antibody

## Preparation

Antigen: Recombinant protein of human Collagen I

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, ICC/IF, FC, IP

## Dilutions


WB: 1/1000

;ICC/IF: 1/50

;FC: 1/50

;IP: 1/20

## Validations



Western blot detection of Collagen I in A549,HL-60,U2OS using Collagen I Rabbit mAb(1:1000 diluted)