

Anti-Phospho-Vimentin (Ser56) Rabbit antibody



Catalog Number: 168705

Product name

Rabbit IgG

Anti-Phospho-Vimentin (Ser56) Rabbit antibody

Applications

Specificity

WB, ICC/IF, FC, IP

Human

Dilutions

Antibody description

WB: 1/1000

Phospho-Vimentin (Ser56) Rabbit monoclonal antibody

ICC/IF: 1/200

FC: 1/20

Preparation

IP: 1/20

Antigen: A synthetic phosphopeptide corresponding to residues surrounding Ser56 of human Vimentin

Validations

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Western blot detection of Phospho-Vimentin (Ser56) in U251MG cell lysates using Phospho-Vimentin (Ser56) Rabbit mAb(1:500 diluted). Predicted band size:54kDa. Observed band size:54kDa.

