

# DRD1 antibody



Catalog Number: 110015

---

## Product name

DRD1 antibody

## Immunogen

[Human DRD1 Recombinant protein \(GST tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

DRD1 Rabbit Polyclonal antibody. Positive WB detected in SH-SY5Y cells. Observed molecular weight by Western-blot: 50 kDa

## Preparation

This antibody was obtained by immunization of DRD1 recombinant protein (Accession Number: NM\_000794). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:500-1:5000

## Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:110015(DRD1 antibody) at dilution of 1:1000