

# DR5-Specific antibody

Catalog Number: 107015

## Product name

DR5-Specific antibody

## Specificity

Human; other species not tested.

## Antibody description

DR5-Specific Rabbit Polyclonal antibody. Positive WB detected in HL-60 cells, PC-3 cells. Observed molecular weight by Western-blot: 58-60 kDa, 45-50 kDa, 32 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_003842). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions


Recommended Dilution:

WB: 1:200-1:1000

## Validations



WB result of Catalog No:110013 (TNFRSF10B, DR5 antibody) with PC-3 cell, untreated or treated



with tunicamysin.

HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:110013(TNFRSF10B-Specific antibody) at dilution of 1:200